ORALMEDIC®

Instructions for use

Barrier Solution

One application of **ORALMEDIC®** stops the pain and discomfort of a mouth ulcer in seconds. A single use prefilled **ORALMEDIC** swab delivers the amount of liquid material required to treat one typical mouth ulcer. It works by forming a protective barrier over the ulcer to provide optimum healing conditions for natural healing.

ORALMEDIC is for the treatment of ulcerations in the tissue of the mouth. It is not intended or effective in the treatment of cold sores, fluid-filled blisters or vesicles in the mouth or on surface of the lips.

WARNING: Do not use if allergic to sulfur.

Instructions:

Important - Please read carefully

1. Dry affected area using the cotton drying swab.

2. Hold the applicator at arms length, away from the eyes, with colored ring up. Gently bend the colored ring to one side until it snaps. **ORALMEDIC** will flow to the opposite end.

3. Apply to the affected area for **5-10 seconds**. Use a rolling motion to completely cover the entire ulcer bed and ulcer rim. A 'stinging' sensation is experienced immediately upon application.

4. **Rinse mouth thoroughly with water and spit out.**

5. Return used swab to package. Dispose of properly.

An effective treatment is evident when a "white patch" is seen in and around the ulcer. This white patch may flake off during the healing period.

Note: If the ulcer pain remains after rinsing with water, it is an indication that some part of the ulcer was not covered. Repeat application immediately, one more time only, following the directions above. Thereafter do not retreat again.

Important Safety information: Do not use if allergic to materials that contain sulfur in any form. Do not use if you are pregnant, or on children under the age of 12. Avoid contact with the eyes. If eye exposure occurs, immediately remove any contact lenses, irrigate eyes for at least 15 minutes with water and seek immediate medical advice.

May cause local irritation or mild swelling when applied. If excess irritation occurs during use, rinse thoroughly with plenty of water. If condition persists, consult your doctor or dentist.

Avoid contact with clothing as it may discolor or damage.

Keep out of reach and sight of children.

Store at room temperature in a dry place.

If used as directed the product cannot be reused. After the swab is activated and liquid is applied, the oral cavity is rinsed with water and spat out, with the used swab properly disposed of.

Ingredients: Hydroxybenzenesulfonic Acid, Hydroxymethoxybenzenesulfonic Acid, Sulfuric Acid, Water, FD&C Red 40.

This pack contains: 3 single use applicators each containing 0.2ml of solution. 3 drying swabs.

EPIEN Medical, Inc. 4225 White Bear Parkway, Ste 600
St. Paul, MN 55110 USA www.epien.com 651-653-3380

ORALMEDIC® is a registered trademark of EPIEN Medical, Inc.

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